## Log file saving in SuperPro 7100

In Superpro 7100, user can save the operations including MAC address read from Gainspan module to a text file (log file).

Click File- Option->Log file in Superpro 7100 menu bar, following window will pop up:

Log File	×
Option Never New Append	OK Cancel
Filename: C:\Users\frys\Desktop\GS.txt	

**Never:** Log file feature is not used and current log information is never saved.

**New:** Create a new log file and overwrite the contents of any such log file with same name.

**Append:** Add information into the current log file. If a lof file already exists, then add new information at the end of current log file.

In following example, we will show how to save current log to GS.txt log file on desktop.

Step 1: Load Project file GS2011.prj or any other project by clicking Project->Load project

**Step 2:** Click File- Option->Log file in Superpro 7100 menu bar, following window will pop up:

Lo	og File	x
Option Never New Append	OK Cancel	
Filename:		

## **Step 3:** Select Append and browse to a file GS.txt on desktop and click OK

Log File	×
Option O Never O New	OK Cancel
Append Filename: C:\Users\frys\Desktop\GS.txt	

**Step4:** Whenever user clicks save (as shown below) or closes SuperPro 7100 software, log information will be saved to GS.txt.

Site #1
Log window
: 00
: 00
: 00
MAC Addr_STK2: 00
: 00
: 00
: 00
: 00
: 00
MAC Addr_STK3: 00
: 00
: 00
: 00
: 00
MAC Addr_STK4: 00
: 00
: 00
: 00
FirmWare_Select:: CURRENT_PREVIOUS_FACTORY VERSION Operation Option
Pins Check: «
Remove interval(s): 4
Insert interval(s): 2
Verify mode: Once with Vcc
Unit Setting: SKT1: / SKT2: / SKT3: / SKT4: ×
Save

For further information, please contact <a href="mailto:techsupport@xeltek.com">techsupport@xeltek.com</a>